

Navina™ Smart product information

Intended use

The Navina Smart is intended for transanal irrigation by instilling water up into the lower part of the colon through a rectal catheter.

Indications for use

Navina Smart is indicated to help adults and children from 3 years who suffer from fecal incontinence, chronic constipation and/or time-consuming bowel management. By instilling water up into the lower part of the colon, the peristaltic muscles in the bowel can be triggered and start to evacuate the content of the lower colon and rectum.

- Always consult a healthcare professional specialized in TAI before using Navina Systems.
- Only use Navina Systems for its intended use, as described in the instruction manual.
- Navina Systems is for a single user and should not be shared with other people.
- The Navina catheter regular is for adult use only.

Contraindications

You should not use Navina Smart system if you have one or more of the following:

- Known anal or colorectal stenosis
- Active inflammatory bowel disease.
- Acute diverticulitis.
- Colorectal cancer.
- Ischemic colitis.
- Within three months of anal or colorectal surgery.
- You are within 4 weeks of previous endoscopic polypectomy.

As the list of contraindications may not be exhaustive, healthcare professionals will always consider individual user factors as well.

Navina Smart Systems

Reference	Type	Specification
69009	Navina Smart System Regular For adult use only	<ul style="list-style-type: none"> 1 Navina smart control unit 1 Power adapter 1 Cable 1 Water container 1 Tube set 2 Rectal catheters Regular 1 Navina case 1 Accessory set
69010	Navina Smart System Small From 3 years	<ul style="list-style-type: none"> 1 Navina Smart control unit 1 Power adapter 1 Cable 1 Water container 1 Tube set 2 Rectal catheters Small 1 Navina case 1 Accessory set
69019	Navina Smart System Cone From 3 years	<ul style="list-style-type: none"> 1 Navina Smart control unit 1 Power adapter 1 Cable 1 Water container 1 Tube set 2 Cones 1 Navina case 1 Accessory set



Navina Smart System Consumables

Reference	Type	Specification	Technical details
69003	Navina Consumable Set Regular	1 Water container 15 Rectal catheters Regular*	15 uses, 1.5 liter, foldable, double scale for visibility. Hydrophilic, single use catheters.
69004	Navina Consumable Set Small	1 Water container 15 Rectal catheters Small	15 uses, 1.5 liter, foldable, double scale for visibility. Hydrophilic, single use catheters.
69016	Navina Consumable Set Cone	1 Water container 15 Cone catheters	15 uses, 1.5 liter, foldable, double scale for visibility. Hydrophilic, single use catheters.
68940	Navina Catheter Set Regular*	10 catheters	Hydrophilic, single use catheters.
68943	Navina Catheter Set Small	10 catheters	Hydrophilic, single use catheters.
69011	Navina Tube Set	1 Catheter/Cone tube 1 Water container tube	100 uses 100 uses
69012	Navina Accessories	2 Grip rings (1 vertical, 1 horizontal) 1 Lanyard 1 Positioning strap 1 Navina Smart positioning clip	
69014	Navina Smart Control Unit	1 Navina Smart control unit 1 Power adaptor 1 Cable	400 uses or two years from the date of the first use
69015	Navina Smart Power Adapter and Cable	1 Power adaptor 1 Cable	none

* For adult use only

All Navina Smart components except for Navina Smart control unit, power adaptor and cable should be disposed of as household waste.

The Navina Smart control unit, power adaptor and cable contain recyclable materials. Do not dispose these product as unsorted municipal waste. Please contact the local municipality for the nearest collection point. For Europe: Separate disposal for electrical and electronic waste according to the Waste Electrical and Electronic Equipment (WEEE Directive).

wellspect.com

