



The funnel has a different color depending on the diameter of the catheter.

Adhesive tab on back for reliable hanging.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>1,2</sup>

URINARY CATHETER  
CATEDE URINARIO

Smooth eyelets.<sup>3,4</sup>



Nelaton Tiemann



User-friendly packaging.

### Product description

LoFric®, single use, hydrophilic urinary catheter, intended for intermittent urinary catheterization.

### Sterilization and shelf life

Sterilized by ethylene oxide gas. Use before expiry date on package.

### Package

30 catheters per customer box. 2 or 4 customer boxes per shipping box.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP) and sodium chloride (NaCl).

#### Catheter:

Polyolefin-based elastomer (POBE).

#### Item packaging:

Lacquered paper, laminated film of Polypropylene (PP)/Polyamide (PA) / Polyethylene (PE).

### Product range

Reference	Type	Length	Size
400 08 XX*	LoFric® Nelaton (male)	40 cm	CH 08
400 10 XX*	LoFric® Nelaton (male)	40 cm	CH 10
400 12 XX*	LoFric® Nelaton (male)	40 cm	CH 12
400 14 XX*	LoFric® Nelaton (male)	40 cm	CH 14
400 16 XX*	LoFric® Nelaton (male)	40 cm	CH 16
400 18 XX*	LoFric® Nelaton (male)	40 cm	CH 18
400 20 XX*	LoFric® Nelaton (male)	40 cm	CH 20
400 22 XX*	LoFric® Nelaton (male)	40 cm	CH 22
401 06 XX*	LoFric® Nelaton (paediatric)	20 cm	CH 06
401 08 XX*	LoFric® Nelaton (paediatric)	20 cm	CH 08
401 10 XX*	LoFric® Nelaton (paediatric)	20 cm	CH 10
402 08 XX*	LoFric® Nelaton (paediatric)	30 cm	CH 08
402 10 XX*	LoFric® Nelaton (paediatric)	30 cm	CH 10
403 08 XX*	LoFric® Nelaton (female)	20 cm	CH 08
403 10 XX*	LoFric® Nelaton (female)	20 cm	CH 10
403 12 XX*	LoFric® Nelaton (female)	20 cm	CH 12
403 14 XX*	LoFric® Nelaton (female)	20 cm	CH 14
403 16 XX*	LoFric® Nelaton (female)	20 cm	CH 16
403 18 XX*	LoFric® Nelaton (female)	20 cm	CH 18
404 08 XX*	LoFric® Nelaton (female)	15 cm	CH 08
404 10 XX*	LoFric® Nelaton (female)	15 cm	CH 10
404 12 XX*	LoFric® Nelaton (female)	15 cm	CH 12
404 14 XX*	LoFric® Nelaton (female)	15 cm	CH 14
405 10 XX*	LoFric® Tiemann (male)	40 cm	CH 10
405 12 XX*	LoFric® Tiemann (male)	40 cm	CH 12
405 14 XX*	LoFric® Tiemann (male)	40 cm	CH 14
405 16 XX*	LoFric® Tiemann (male)	40 cm	CH 16
405 18 XX*	LoFric® Tiemann (male)	40 cm	CH 18

\*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

## References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

### **Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

#### LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

